

# Cell Depletion

with Akadeum Microbubbles

## 1. Prepare your cells

- 1.1 Homogenize tissue and lyse RBCs as needed.
- 1.2 Centrifuge cell suspension (5 min, RT, 400 × g), aspirate supernatant, and wash once with Separation Buffer.
- 1.3 Centrifuge samples (5 min, RT, 400 × g) and aspirate supernatant.
- 1.4 Resuspend cells in Separation Buffer and transfer to 5 mL Eppendorf tubes.

### Recommended supplies:

1. Cold (2-8 °C) Separation Buffer (Ca<sup>2+</sup> and Mg<sup>2+</sup> free PBS containing 2 mM EDTA and 0.5% biotin free BSA)
2. Low retention 1 mL pipet tips (VWR Part #89174-530)
3. Centrifuge with swinging bucket rotor

## 2. Label cells

- 2.1 Add biotin-antibody cocktail. Mix briefly.
- 2.2 Incubate sample for 20 min at 4°C.
- 2.3 Add Separation Buffer and centrifuge (5 min, RT, 400 × g). Aspirate supernatant.
- 2.4 Resuspend cells in Separation Buffer.

## 3. Bind microbubbles

- 3.1 Prepare microbubbles by resuspending in solution (microbubble mixture should be a homogenous white solution, *i.e.* look like milk). Vigorously mix or pipet. Immediately proceed to next step.
- 3.2 Add microbubbles to first sample.
- 3.3 Set the pipet volume to ~70% of the total sample volume (cell suspension + microbubbles) and mix microbubbles with gentle trituration for 30 pipet strokes using low retention 1000 µL pipet tip.  
*Note: This pipet setting ensures adequate mixing of microbubbles and cells for binding.*
- 3.4 Immediately add Separation Buffer.  
*Note: This facilitates separation of microbubbles from the cell pellet.*
- 3.5 Repeat steps 3.1–3.4 for remaining samples making sure to **resuspend microbubbles before each sample.**

## 4. Separate cells

- 4.1 Centrifuge samples (5 min, RT, 400 × g).  
*Note: Use of a swinging bucket rotor for this step facilitates microbubble aspiration.*
- 4.2 Aspirate off white microbubble layer and supernatant. Take care not to aspirate cell pellet.
- 4.3 Resuspend cells in Separation Buffer. Transfer cells to a new tube if desired due to residual microbubbles stuck to tube.
- 4.4 Continue to flow cytometry or other downstream processing or handling.