

# Mouse Naïve CD4+ T Cell Isolation Kit (Cat.12210-150)

## BACS™ Microbubbles Protocol

Isolate untouched Naïve CD4+ T cells from mouse splenocytes via Buoyancy Activated Cell Sorting (BACS). This kit can be used to target and remove non-naïve CD4+ T cells with antibodies recognizing CD8a, CD11b, CD11c, CD19, CD24, CD25, CD44, CD49b, CD105, Gr1, Ter119, and TCR $\gamma\delta$ . Isolated naïve CD4+ T cells are suitable for flow cytometry, molecular assays, activation and expansion, cell culture, and other functional studies. Processing capacity  $1 \times 10^9$  cells.

| Name   | Format   | Quantity     | Storage |
|--|--|--------------|---------|
| BACS™ Streptavidin Microbubbles                  | In buffer with 0.09% sodium azide.   | 10.5 mL      | 2-8 °C  |
| Mouse Naïve CD4+ T Cell Biotin Antibody Cocktail | Monoclonal antibodies in PBS with sodium azide.  | 1050 $\mu$ L | 2-8 °C  |
| Separation Buffer                                | Ca <sup>2+</sup> and Mg <sup>2+</sup> –free PBS containing 2 mM EDTA and 0.5% biotin-free BSA. | 200 mL       | 2-8 °C  |
| 5 mL Tubes                                       | Bag of tubes   | 20 tubes     | RT      |

### Additional Supplies:



- 1 20 rpm tube rotator for mixing (e.g., Thermo Scientific cat#: 88881002)
- 2 Centrifuge (swinging bucket rotor strongly preferred)
- 3 Vacuum aspirator
- 4 30  $\mu$ m cell strainer (optional)

### Before You Begin:

- ▶ This protocol has been optimized for splenocytes as the starting material. If working with other sample types, please contact techsupport@akadeum.com.
- ▶ For optimal results, homogenize mouse spleens and lyse red blood cells in the sample prior to separation.
- ▶ Separation Buffer is azide-free. Cell isolation should be conducted under aseptic conditions.
- ▶ For optimal results, prior to cell separation, filter samples through a 30  $\mu$ m cell strainer to obtain a single-cell suspension.
- ▶ For tips on how to vacuum aspirate the BACS™ Microbubble layer, see video: <https://www.akadeum.com/videos/aspiration>
- ▶ This protocol is designed for starting samples containing  $1 \times 10^7$  –  $24 \times 10^7$  total cells. Samples with  $> 24 \times 10^7$  should be divided across multiple tubes. For samples  $< 1 \times 10^7$ , please contact techsupport@akadeum.com

### Representative Naïve CD4+ T Cell Isolation:

Untouched naïve T cells were isolated from splenocytes. Isolated cells were labeled with CD62L-PerCPCy5.5 and CD4-APCFire750. The fluorescently labeled cells were analyzed by flow cytometry. Debris and dead cells were excluded from analysis.



## Experimental Setup:

| Sample Size               | Tube Size | Sample Volume (Step 2)        | Antibody Cocktail (Step 4)    | BACSTM Microbubbles (Step 6)  | Final Volume (Step 7) |
|---------------------------|-----------|-------------------------------|-------------------------------|-------------------------------|-----------------------|
| (1x10 <sup>7</sup> cells) |           | per (1x10 <sup>7</sup> cells) | per (1x10 <sup>7</sup> cells) | per (1x10 <sup>7</sup> cells) | Separation Buffer     |
| 1 - 7                     | 1.5 mL    | 30 µL                         | 10 µL                         | 100 µL                        | Fill to 1.2 mL        |
| > 7 - 24                  | 5.0 mL    | 30 µL                         | 10 µL                         | 100 µL                        | Fill to 4.0 mL        |

## Prepare Cells:

- Count and wash cells.
- Resuspend cell pellet in 30 µL of Separation Buffer per 1 x 10<sup>7</sup> cells, as indicated in the table above.
- Transfer cell suspension to a 1.5 or 5 mL tube, as indicated in the table above. Divide or aliquot sample to be within the cell number ranges indicated in the table above.
- Add 100 µL of BACSTM Microbubbles per 1 x 10<sup>7</sup> total cells to the labeled sample as indicated in the table above.
- Add Separation Buffer to achieve a final volume of 1.2 or 4.0 mL, as indicated in the table above.
- Mix samples on a rotator at 20 rpm for 10 min at room temperature (or at 4°C).

## Label Cells:

- Add 10 µL of Mouse Naïve T CD4+ Cell Isolation Biotin Antibody Cocktail per 1 x 10<sup>7</sup> total cells as indicated in the table above. Gently mix samples and incubate for 10 min at room temperature (or at 4°C).

## Bind BACSTM Microbubbles:

- Resuspend BACSTM Microbubbles by pipetting or inverting by hand.

*Note: It is critical that BACSTM Microbubbles are thoroughly resuspended immediately prior to addition to each sample. Resuspension can be achieved by pipetting with a 1 mL pipette 2-3 times, followed by inverting multiple times to create a homogeneous suspension.*

## Separate Cells:

- Centrifuge samples at 400 x g for 5 min.  
*Note: A swinging bucket rotor centrifuge is recommended.*
- Vacuum aspirate the BACSTM Microbubble layer and supernatant, taking care not to disturb the cell pellet. Once BACSTM Microbubbles have been aspirated, the supernatant may be removed by pipette.  
*Note: For tips on how to remove BACSTM microbubbles, see video: <https://www.akadeum.com/videos/aspiration>*
- Resuspend cell pellet in desired buffer or media and transfer to clean tube.



**AKADEUM**  
LIFE SCIENCES

## Safety Information

For research use only. Not intended for any animal or human therapeutic or diagnostic use. For information regarding hazards and safe handling practices, please consult the Safety Data Sheet.

Visit us here for how-to videos, additional product information, and tech support

